



# C-Cup Multi-Drug Urine Test Cup

Catalogue No. See Box Label

The SAFElife® C-Cup Multi-Drug Urine Test Cup contains competitive binding, lateral flow immunochemotographic assays for qualitative and simultaneous detection of 6-Monocetyl morphine, Amphetamine, Sebacorbital, Buprenorphine, Oxycodone, Cocaine, Caffeine, 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP), Ethyl Glucuronide, Fentanyl, Gabapentin, Hydromorphone, Synthetic Cannabinoids, Ketamine, Kratom, Lysergic acid diethylamide, Methylenedioxymethamphetamine, Methamphetamine, Morphine, Methadone, Methqualone, Opiate, Oxycodone, Phencyclidine, Pregabalin, Propoxyphene, Nortriptyline, Cannabinoids, Tramadol and Alcohol in human urine with below cutoff concentrations and approximate detection time:

Drug (Identifier)	Calibrator	Cut-off Level	Minimum Detection Time	Maximum Detection Time
6-Monocetyl morphine (6-MAM)	6-Monocetyl morphine	10 ng/mL	2 hours	8 hours
Amphetamine (AMP300)	d-Amphetamine	300 ng/mL	2-7 hours	1-2 days
Amphetamine (AMP500)	d-Amphetamine	500 ng/mL	2-7 hours	1-2 days
Sebacorbital (BAR)	d-Amphetamine	1000 ng/mL	2-7 hours	1-2 days
Buprenorphine (BUP)	Buprenorphine	5 ng/mL	2-4 hours	1-3 days
Buprenorphine (BUP10)	Buprenorphine	10 ng/mL	4 hours	1-3 days
Oxazepam (BZO100)	Oxazepam	100 ng/mL	2-7 hours	1-2 days
Oxazepam (BZO200)	Oxazepam	200 ng/mL	2-7 hours	1-2 days
Oxazepam (BZO300)	Oxazepam	300 ng/mL	2-7 hours	1-2 days
Cocaine (COC100)	Benzoylecgonine	100 ng/mL	1-4 hours	2-4 days
Cocaine (COC200)	Benzoylecgonine	150 ng/mL	1-4 hours	2-4 days
Cocaine (COC300)	Benzoylecgonine	300 ng/mL	1-4 hours	2-4 days
Cotinine (COT)	Cotinine	200 ng/mL	2-8 hours	1-7 days
EDDP100	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	100 ng/mL	3-8 hours	1-3 days
EDDP300	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300 ng/mL	3-8 hours	1-3 days
Ethyl (EG300)	Glucuronide	300 ng/mL	1-2 hours	Up to 3+ days
Ethyl (EG500)	Glucuronide	500 ng/mL	1-2 hours	Up to 3+ days
Fentanyl (FTY)	Nortriptyline	20 ng/mL	1-4 hours	1-3 days
Gabapentin (GAB)	Gabapentin	2000 ng/mL	5-7 hours	Up to 2 days
Hydromorphone (HMO)	Hydromorphone	300 ng/mL	4-6 hours	1-2 days
Synthetic Cannabinoid (K2)	JWH-073 Butanoic Acid	50 ng/mL	8-12 hours	Up to 5+ days
Ketamine (KET300)	Ketamine	300 ng/mL	2-4 hours	2-3 days
Ketamine (KET1000)	Ketamine	1000 ng/mL	2-4 hours	2-3 days
Kratom (KRA100)	Mitragynine	100 ng/mL	7 hours	5-6 days
Kratom (KRA300)	Mitragynine	300 ng/mL	7 hours	5-6 days
Lysergic acid diethylamide (LSD)	Lysergic acid diethylamide	20 ng/mL	2.5 hours	Up to 5+ days
Methylenedioxymethamphetamine (MDMA)	3,4-Methylenedioxymethamphetamine (MDMA)	500 ng/mL	2-7 hours	2-4 days
Methamphetamine (MET300/mAMP300)	Di(-)-Methamphetamine	300 ng/mL	2-7 hours	2-4 days
Methamphetamine (MET500/mAMP500)	Di(-)-Methamphetamine	500 ng/mL	2-7 hours	2-4 days
Methamphetamine (MET1000/mAMP1000)	Di(-)-Methamphetamine	1000 ng/mL	2-7 hours	2-4 days
Morphine (MOP100/OP100)	Morphine	100 ng/mL	2 hours	2-3 days
Morphine (MOP300/OP300)	Morphine	300 ng/mL	2 hours	2-3 days
Methadone (MTD200)	Methadone	200 ng/mL	3-8 hours	1-3 days
Methadone (MTD300)	Methadone	300 ng/mL	3-8 hours	1-3 days
Methqualone (MQ)	Methqualone	300 ng/mL	6-8 hours	Up to 7+ days
Opiate (OPI2000)	Morphine	2000 ng/mL	2 hours	2-3 days
Oxycodone (OXY)	Oxycodone	100 ng/mL	4 hours	1-3 days
Phencyclidine (PCP)	Phencyclidine	25 ng/mL	4-6 hours	7-14 days
Pregabalin (PGB)	Pregabalin	500 ng/mL	6-8 hours	1-3 days
Propoxyphene (PPX)	Propoxyphene	300 ng/mL	2 hours	2-3 days
Nortriptyline (TCA)	Nortriptyline	1000 ng/mL	8-12 hours	2-7 days
Cannabinoids (THCS)	11-nor-Δ <sup>9</sup> -THC-9-COOH	15 ng/mL	2 hours	Up to 5+ days
Cannabinoids (THCS2)	11-nor-Δ <sup>9</sup> -THC-9-COOH	25 ng/mL	2 hours	Up to 5+ days
Cannabinoids (THC40)	11-nor-Δ <sup>9</sup> -THC-9-COOH	40 ng/mL	2 hours	Up to 5+ days
Cannabinoids (THC50)	11-nor-Δ <sup>9</sup> -THC-9-COOH	50 ng/mL	2 hours	Up to 5+ days
Tramadol (TRA 100)	Tramadol	100 ng/mL	8-12 hours	3-7 days
Tramadol (TRA 200)	Tramadol	200 ng/mL	8-12 hours	3-7 days
Tramadol (TRA 1000)	Tramadol	1000 ng/mL	8-12 hours	3-7 days
Alcohol (ETOH)	Alcohol	0.04 g/dL	-	-

Configurations of the SAFElife® C-Cup Multi-Drug Urine Test Cup can consist of any combination of the above listed drug analytes.

It is intended for forensic use only.

It is not intended to distinguish between prescription use or abuse of these drugs. Professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS) is the recommended confirmatory method.

## WARNINGS AND PRECAUTIONS

- The test kit is for external use only.
- Discard after first use. The test kit cannot be used more than once.
- Do not use the test kit beyond expiration date.
- Do not use the test kit if the pouch is punctured or not well sealed.
- Keep out of the reach of children.

## Approximate Alcohol Concentration

0.0 mg/100 mL	40 mg/100 mL	80 mg/100 mL	200 mg/100 mL
	(0.04%)	(0.08%)	(0.2%)

## Invalid

The test should be considered invalid if only the edge of the reaction pad turned color that might be ascribed to insufficient sampling. Another test should be run to re-evaluate the specimen. If test still fails, please contact the distributor, with the lot number.

## What Is the False Positive Test?

The definition of the false positive test would be an instance where a substance is identified incorrectly by the SAFElife® C-Cup Multi-Drug Urine Test Cup. The most common causes of the false positive test are cross reactants. Certain foods and medicines, diet plan drugs and nutritional supplements may cause the false positive test result.

## What Is the False Negative Test?

The definition of the false negative test is that the initial drug is present but isn't detected by the SAFElife® C-Cup Multi-Drug Urine Test Cup. If the specimen is diluted or adulterated, it may cause false negative result.

If suspect someone is taking drugs but get the negative test results, please test again at another time.

## TEST LIMITATIONS

- This test kit has been developed for testing urine specimen only. No other fluids have been evaluated. DO NOT use it to test anything other than urine.
- Adulterated urine specimen may produce false results. Strong oxidizing agents such as bleach (hypochlorite) can oxidize drug analytes. If a specimen is suspected of being adulterated, obtain a new specimen.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause false results.
- This test is a qualitative screening assay. It is not designed to determine the quantitative concentration of drugs or the level of intoxication.

## SUMMARY

**6-Monocetyl morphine (6-MAM)**  
Heroin is rapidly metabolized in the body. The half-life in blood is only 3-19 minutes. It is degraded by esterase in the body to 6-monocetyl morphine (hereinafter abbreviated as 6-MAM) and the molecular formula is C<sub>21</sub>H<sub>27</sub>NO<sub>5</sub>. 6-MAM is detected in the body to form morphine, and morphine cannot be acetylated to form 6-MAM in vivo. 6-MAM in the human body is only derived from the metabolism of heroin. Therefore, the U.S. Department of Health and Human Services (DHHS) recommended 6-monocetyl morphine as a specific test for heroin abuse.

**Amphetamine (AMP)**  
Amphetamine and the structurally related "designer" drugs are sympathomimetic amines whose biological effects include potent central nervous system (CNS) stimulation, anorectic, hyperthermic, and cardiovascular properties. They are usually taken orally, intravenously, or by smoking. Amphetamine is readily absorbed from the gastrointestinal tract and are then either deactivated by the liver or excreted unchanged in the urine with a half-life of about 12 hours. It can be detected in the urine for 1 to 2 days after use. Amphetamine is metabolized to deaminated (hippuric and benzoic acids) and hydroxylated metabolites. Methamphetamine is partially metabolized to amphetamine and its major active metabolite. Amphetamines increase the heart rate and blood pressure, and suppress the appetite. Studies indicate that heavy abuse may result in permanent damage to certain essential nerve structural in the brain.

**Sebacorbital (BAR)**  
Sebacorbital are a class of central nervous system depressants. They have a wide range of half-life of 2 to 40 hours and can be detected in the urine for 1 to 4 days after use. Phenobarbital is a long acting barbiturate derivative that has been used as a daytime sedative and very extensively as an anticonvulsant. Phenobarbital and sebacorbital are excreted in the urine unchanged. Abuse of barbiturates can lead not only to impaired motor coordination and mental disorder, but also to respiratory collapse, coma and even death. Barbiturates are taken orally, rectally, or by intravenous and intramuscular injections. Short-acting barbiturates will generally be excreted in urine as metabolites, while the long-acting barbiturates will primarily appear unchanged.

**Buprenorphine (BUP)**  
Buprenorphine is a partial analgesic often used in the treatment of opioid addiction. The drug is sold under the trade names Subutex®, Buprenex®, Temgesic® and Suboxone®; all of which contain Buprenorphine HCl alone or in combination with Naloxone HCl. Buprenorphine is used as a substitution treatment for opioid addicts. A substitution treatment is a form of medical care offered to opiate addicts (primarily heroin addicts) based on a similar or identical substance to the drug normally used. In substitution therapy, Buprenorphine is as effective as Methadone but demonstrates a lower level of physical dependence. The plasma half-life of Buprenorphine is 2-4 hours. While complete elimination of a single dose of the drug can take as long as 6-9 days, the detection window for the parent drug in urine is thought to be approximately 3 days.

**Oxazepam (BZO)**  
Benzodiazepines are the most widely used anxiolytic drugs. They are used extensively as anti-anxiety agents, hypnotics, muscle relaxants and anti-convulsants. They are taken orally or sometimes by injection and have a wide range of half-life from 2 to 40 hours. They can generally be detected for 1 to 2 days after Benzodiazepines use. Benzodiazepines are metabolized in the liver. Some Benzodiazepines and their metabolites are excreted in the urine. Their use can result in drowsiness and/or confusion. Benzodiazepines potentiate alcohol and other CNS depressants. Psychological and physical dependence on benzodiazepines can develop if high doses of the drug are given over a prolonged period.

**Cocaine (COC)**  
Cocaine derived from leaves of coca plant, is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using cocaine are euphoria, confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Cocaine is excreted in urine primarily as benzoylecgonine in a short period of time.

**Morphine (MOP/OP)**  
The opiates such as heroin, morphine, and cocaine are derived from the resin of opium poppy. The principal metabolites of opiates are morphine, morphine-3-glucuronide morphine and cocaine with a half-life of about 3 hours. Heroin is quickly metabolized to morphine. Thus, morphine and morphine glucuronide might be detected in the urine. Their use can result in drowsiness and/or confusion. Benzodiazepines potentiate alcohol and other CNS depressants. Psychological and physical dependence on benzodiazepines can develop if high doses of the drug are given over a prolonged period.

**Methamphetamine (MET/mAMP)**  
Methamphetamine is a potent sympathomimetic agent with therapeutic applications. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and a sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behavior, and cardiac dysrhythmia. The pattern of psychomotor activity may appear at half-life of about 15 hours and is excreted in urine as amphetamine and oxidized as deaminated and hydroxylated derivatives. However, 40% of methamphetamine is excreted unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

**Cocaine (COC)**  
Cocaine derived from leaves of coca plant, is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using cocaine are euphoria, confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Cocaine is excreted in urine primarily as benzoylecgonine in a short period of time.

**Cotinine (COT)**  
Cotinine is an alkaloid found in tobacco and is also a major metabolite of Nicotine, which produces stimulation of the autonomic ganglia and central nervous system when in humans. Nicotine is found in tobacco products such as cigarettes, tobacco chew, and nicotine patches or gums. It is an addictive substance and is poisonous in a large amount. In addition to addiction, some of the other substances within tobacco products, such as carbon monoxide or tar, are dangerous to the body and can lead to medical conditions such as emphysema, lung cancer, and heart disease. In a 24-hour urine, approximately 5% of a nicotine dose is excreted as unchanged drug with 10% as cotinine and 35% as hydroxycotinine. The concentrations of other metabolites are believed to account for less than 5%. While Cotinine is thought to be an inactive metabolite, its elimination profile is more stable than that of Nicotine which is largely urine PH dependent. Cotinine is stable in body fluids and has a relatively long half-life of approximately 17 hours, and is typically detectable for several days after the use of tobacco. Therefore, the presence of cotinine in the detection of Cotinine is less dependent on the time of sampling than that of Nicotine. Nicotine and Cotinine are rapidly eliminated by the kidney; the window of detection for cotinine in urine at a cutoff level of 200

ng/mL is expected to be up to 2-3 days after nicotine use.

## EDDP

EDDP (2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine) is the primary metabolite of methadone. Methadone is a synthetic analgesic drug that is originally used in the treatment of narcotic addicts. The detection of EDDP is more beneficial than traditional methadone screening since EDDP exists only in urine fluids, indicating that ingested methadone. The tolerance of specimens by spiking the urine with methadone can be prevented. Secondly, renal clearance of EDDP is not affected by urinary pH, therefore the EDDP test provides a more accurate result of methadone ingestion than the methadone parent screening.

## Ethyl Glucuronide (EtG)

Ethyl Glucuronide is a direct metabolite of alcohol. Presence in urine may be used to detect recent alcohol intake, even after alcohol is no longer measurable. Traditional laboratory methods detect the actual alcohol in the body, which reflects current intake within the past few hours (depending on how much was consumed). The presence of EtG in urine is a definitive indicator that it can be detected in the urine for 2 to 4 days after drinking alcohol, even alcohol is eliminated from the body. Therefore, EtG is a more accurate indicator of the recent intake of alcohol than measuring for the presence of alcohol itself. The EtG test can aid in the diagnosis of drunk driving and alcoholism, which has important significance in the forensic identification and medical examination.

## Fentanyl (FTY)

Fentanyl is a potent, synthetic narcotic analgesic with a rapid onset and short duration of action. It was first synthesized by Janssen Pharmaceutische (Belgium) in the late 1950s, and it is approximately 100 times more potent than morphine. Fentanyl is a strong agonist at the μ-opioid receptor. Historically it has been used to treat breakthrough pain and is commonly used in pre-procedures as a pain reliever as well as an anesthetic in combination with a benzodiazepine. Fentanyl is frequently given intrathecally as part of spinal anesthesia or epidurally for epidural anesthesia and analgesia.

## Gabapentin (GAB)

Gabapentin (GAB), sold under the brand name Neurontin, is a medication used to treat epilepsy, neuropathic pain, hot flashes, and restless legs syndrome. In epilepsy, it may be used for those with partial onset (approximately 50% of the recommended dose) or for the treatment of tonic-clonic seizures. It is used in diabetic neuropathy, postherpetic neuralgia, and central neuropathic pain. It is also used to relieve nerve pain following shingles (a painful rash due to herpes zoster infection) in adults. The most common side effects of gabapentin include dizziness, fatigue, weakness, ataxia, peripheral edema (swelling of extremities), myalgias, and tremor. Serious side effects may include an increased risk of suicide, aggressive behavior, and drug reaction with eosinophilia and systemic symptoms.

## Hydromorphone (HMO)

Hydromorphone, also known as dihydromorphone or dihydromorphine, is a semi-synthetic strong narcotic analgesic. Its chemical structure is similar to morphine, but its coordination, loss of coordination, trance-like ecstatic states, a sense of euphoria and visual distortions. It has many street names, such as "angel dust" and "crystal cyclone", etc. pencyclidine can be administered orally, by nasal ingestion, smoking, or by intravenous injection. It is metabolized in the liver and excreted through the kidneys in urine with a half-life of about 12 hours. Suboxone is a combination of buprenorphine and naloxone.

## Phencyclidine (PCP)

Phencyclidine is an arylcyclohexylamine that was originally used as an anesthetic agent and a veterinary tranquilizer. Phencyclidine can produce hallucinations, lethargy, disorientation, loss of coordination, trance-like ecstatic states, a sense of euphoria and visual distortions. It has many street names, such as "angel dust" and "crystal cyclone", etc. pencyclidine can be administered orally, by nasal ingestion, smoking, or by intravenous injection. It is metabolized in the liver and excreted through the kidneys in urine with a half-life of about 12 hours. Suboxone is a combination of buprenorphine and naloxone.

## Pregabalin (PGB)

Pregabalin (PGB) is a synthetic analogue of γ-aminobutyric acid (GABA), which is similar to gabapentin in structure and action, and has antiepileptic, analgesic and anxiolytic activities. About 98% of pregabalin is recovered in the urine as the active drug after radiation labeling. Therefore, pregabalin abuse can be determined directly by measuring the amount of pregabalin in urine.

## Kratom (KRA)

Kratom (Kratom) (Mitragyna speciosa) is a plant indigenous to Thailand and Southeast Asia. Kratom leaves produce complex stimulant and opioid-like analgesic effects. In Asia, it is often used to stave off fatigue and to manage pain, diarrhea, cough, and opioid withdrawal. Recently, kratom has become widely available in the United States and Europe by means of smoke shops and the Internet. The clinical manifestations of kratom are not well defined and studies are limited, but its safety profile has become a national and international concern, primarily due to excessive consumption being linked to increase in hospital visits for septic injury, seizures, coma, and death. The main active ingredients include Mitragynine and 7-Hydroxymitragynine, which can be detected in urine up to 72 hrs (1-3).

## Lysergic acid diethylamide (LSD)

Lysergic acid diethylamide (LSD) is a white powder or colorless liquid that is a strong semi-artificial hallucinogen. LSD is manufactured from Lysergic acid which occurs naturally in the ergot fungus that grows on wheat and rye. It is a schedule I controlled substance, available in liquid, powder, tablet (microdots), and capsule form. LSD is a non-selective 5-HT<sub>2A</sub> agonist, can exert its hallucinogenic effect by interacting with 5-HT<sub>2A</sub> receptors as a partial agonist and modulating the 5-HT<sub>2A</sub> receptors, producing a marked slowing of the firing rate of serotonergic neurons. LSD can cause the user's senses, feelings, memory, and self-awareness to intensify and change for 6 to 12 hours. In addition to causing mental confusion, LSD can also cause physical pain, with symptoms in the nervous system, cardiovascular, and digestive systems. Most LSD users use marijuana, heroin, or other drugs together.

**Methylenedioxymethamphetamine (MDMA)**  
Methylenedioxymethamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlander, 1990). The most pervasive effect of MDMA, occurring in virtually all people who took a reasonable dose of the drug, was to produce a clenching of the jaws.

**Methamphetamine (MET/mAMP)**  
Methamphetamine is a potent sympathomimetic agent with therapeutic applications. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and a sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behavior, and cardiac dysrhythmia. The pattern of psychomotor activity may appear at half-life of about 15 hours and is excreted in urine as amphetamine and oxidized as deaminated and hydroxylated derivatives. However, 40% of methamphetamine is excreted unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

**Cocaine (COC)**  
Cocaine derived from leaves of coca plant, is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using cocaine are euphoria, confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Cocaine is excreted in urine primarily as benzoylecgonine in a short period of time.

**Cotinine (COT)**  
Cotinine is an alkaloid found in tobacco and is also a major metabolite of Nicotine, which produces stimulation of the autonomic ganglia and central nervous system when in humans. Nicotine is found in tobacco products such as cigarettes, tobacco chew, and nicotine patches or gums. It is an addictive substance and is poisonous in a large amount. In addition to addiction, some of the other substances within tobacco products, such as carbon monoxide or tar, are dangerous to the body and can lead to medical conditions such as emphysema, lung cancer, and heart disease. In a 24-hour urine, approximately 5% of a nicotine dose is excreted as unchanged drug with 10% as cotinine and 35% as hydroxycotinine. The concentrations of other metabolites are believed to account for less than 5%. While Cotinine is thought to be an inactive metabolite, its elimination profile is more stable than that of Nicotine which is largely urine PH dependent. Cotinine is stable in body fluids and has a relatively long half-life of approximately 17 hours, and is typically detectable for several days after the use of tobacco. Therefore, the presence of cotinine in the detection of Cotinine is less dependent on the time of sampling than that of Nicotine. Nicotine and Cotinine are rapidly eliminated by the kidney; the window of detection for cotinine in urine at a cutoff level of 200

ng/mL is expected to be up to 2-3 days after nicotine use.

**Methamphetamine (MET/mAMP)**  
Methamphetamine is a potent sympathomimetic agent with therapeutic applications. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and a sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behavior, and cardiac dysrhythmia. The pattern of psychomotor activity may appear at half-life of about 15 hours and is excreted in urine as amphetamine and oxidized as deaminated and hydroxylated derivatives. However, 40% of methamphetamine is excreted unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

**Cocaine (COC)**  
Cocaine derived from leaves of coca plant, is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using methadone are analgesia, sedation and respiratory depression. Overdose of methadone may cause coma or even death. It is administered orally or intravenously and is metabolized in the liver and excreted in the urine as methadone, EDDP, EMMA and methadone. The kidneys are a major route of methadone excretion. Methadone has a biological half-life of 15 to 60 hours.

**Methadone (MTD)**  
Methadone is a synthetic analgesic drug that is originally used in the treatment of narcotic addicts. Among the psychological effects induced by using methadone are analgesia, sedation and respiratory depression. Overdose of methadone may cause coma or even death. It is administered orally or intravenously and is metabolized in the liver and excreted in the urine as methadone, EDDP, EMMA and methadone. The kidneys are a major route of methadone excretion. Methadone has a biological half-life of 15 to 60 hours.

## Methqualone (MQ)

Methqualone is a sedative that falls outside the benzodiazepine and barbiturate classes. It was once a popular pharmacological agent for sedation and recreational use, but its current use is largely relegated to Africa, particularly South Africa. Because it faced fewer restrictions when it first entered the market, the drug was widely prescribed and perceived as uniquely safe. We now know methqualone can be used recreationally and can cause physical dependence. A lot of lore exists around the drug. In reality, it's not a massively unique drug, and its medical use is limited. The tolerance of specimens by spiking the urine with methadone can be prevented. Secondly, renal clearance of EDDP is not affected by urinary pH, therefore the EDDP test provides a more accurate result of methadone ingestion than the methadone parent screening.

## Opiate (OPI)

Opiate refers to any drug that is derived from the opium poppy, including the natural products, morphine and codeine, and the semi-synthetic drugs such as heroin. Opioid is more general, referring to any drug that acts on the opioid receptor. Opioid analgesics comprise a large group of substances which control pain by depressing the central nervous system. Large doses of morphine can produce higher tolerance levels, physiological dependency in users, and may lead to substance abuse. Morphine is excreted in the urine as morphine-3-glucuronide and is also the major metabolic product of codeine and heroin. Morphine is detectable in the urine for several days after an opiate dose. The test for Morphine 2000 (OPI) of the SAFElife® C-Cup Multi-Drug Urine Test Cup yields a positive result when the morphine in urine exceeds 2000 ng/mL.

## Oxycodone (OXY)

Oxycodone is known as Oxycotin and Roxicodone. It is an ingredient of Percodan, Percocet, Roxicet and Tylox. Oxycodone is a semi-synthetic opiates derived from opium. Like other opiates, Oxycodone is characterized by its analgesic properties, and the tendency for users to form a physical dependency and develop tolerance with extended use. Oxycodone is usually administered in combination with non-opioid analgesics such as acetaminophen and salicylates for the relief of moderate to severe pain. Oxycodone is a central nervous system depressant that may cause drowsiness, dizziness, lethargy, weakness and confusion. Toxicity in an overdose of Oxycodone can lead to stupor, coma, muscle flaccidity, severe respiratory depression, hypotension, and cardiac arrest. Oxycodone is metabolized by H<sub>2</sub> and O<sub>2</sub> demethylation. One of the metabolites, oxycodone, is a potent narcotic analgesic, while the other, noroxycodone, is relatively inactive. From 33 to 61% of a single dose of Oxycodone is excreted in a 24-hour urine collection and consists of 13-19% free Oxycodone, 7-29% glucuronide conjugated Oxycodone, 4-14% glucuronide conjugated oxycodone and an unknown amount of noroxycodone. The detection time window of Oxycodone is 1-3 days following use.

## Phencyclidine (PCP)

Phencyclidine is an arylcyclohexylamine that was originally used as an anesthetic agent and a veterinary tranquilizer. Phencyclidine can produce hallucinations, lethargy, disorientation, loss of coordination, trance-like ecstatic states, a sense of euphoria and visual distortions. It has many street names, such as "angel dust" and "crystal cyclone", etc. pencyclidine can be administered orally, by nasal ingestion, smoking, or by intravenous injection. It is metabolized in the liver and excreted through the kidneys in urine with a half-life of about 12 hours. Suboxone is a combination of buprenorphine and naloxone.

## Pregabalin (PGB)

Pregabalin (PGB) is a synthetic analogue of γ-aminobutyric acid (GABA), which is similar to gabapentin in structure and action, and has antiepileptic, analgesic and anxiolytic activities. About 98% of pregabalin is recovered in the urine as the active drug after radiation labeling. Therefore, pregabalin abuse can be determined directly by measuring the amount of pregabalin in urine.

## Propoxyphene (PPX)

Propoxyphene, a synthetic opiate agonist, is structurally similar to methadone. Propoxyphene is a narcotic analgesic used to relieve mild to moderate pain. The principal metabolites are norendropropoxyphene. The combination usage of propoxyphene, aspirin, acetaminophen or other sedatives can lead cooperative interaction. Abuse of propoxyphene can lead nausea, vomiting, constipation, ileus, heart failure, poisoning, lung dropsy and even death. Propoxyphene is metabolized in the liver and excreted in urine as norendropropoxyphene. Thus the presence of the propoxyphene or its metabolites in the urine indicates propoxyphene use.

## Nortriptyline (TCA)

TCA (Tricyclic Antidepressants) are commonly used for the treatment of depressive disorders. TCA overdoses can result in profound central nervous system depression, cardiotoxicity and anticholinergic effects. TCA overdose is the most common cause of death from prescription drugs. TCAs are taken orally sometimes by injection. TCAs are metabolized in the liver. Both TCAs and their metabolites are excreted in urine mostly in the form of metabolites for up to ten days.

## Cannabinoids (THC)

Cannabinoids are hallucinogenic agents derived from the flowering portion of the hemp plant. The active ingredients in Cannabinoids, THC & Cannabinol can be metabolized and excreted as 11-nor-Δ<sup>9</sup>-tetrahydrocannabinol-9-carboxylic acid with a half-life of 24 hours. They can be detected for 1 to 5 days after use. Smoking is the primary method of use for Cannabinoids/cannabis. Higher doses used by abusers produce central nervous system effects, altered mood and sensory perceptions, loss of coordination, impaired short-term memory, anxiety, paranoia, depression, confusion, hallucinations and increased heart rate. A tolerance to the cardiac and parasympathetic effects can occur, and withdrawal syndrome produces restlessness, insomnia, anorexia and nausea.

## Tramadol (TRA)

Tramadol (2-[(dimethylamino)methyl]-1-[3-(methoxyphenyl) cyclohexan]ol) is used similarly to codeine, to treat moderate to moderate/severe pain. It is a synthetic analog of the phenanthrene alkaloid codeine and, as such, is an opioid and also a pro-drug (codeine is metabolized to morphine, tramadol is converted to O-desmethylen tramadol). Tramadol is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the



THC (50)	Viewer C	+	0	0	2	20	20	100% (91.2% - 100%)
	Viewer B	-	10	13	15	0	0	95% (83.5% - 98.6%)
	Viewer A	+	0	0	1	18	22	100% (84.5% - 100%)
TRA (100)	Viewer C	+	10	12	17	0	0	97.5% (82% - 100%)
	Viewer B	-	0	0	1	18	22	100% (84.5% - 100%)
	Viewer A	+	0	0	1	19	20	97.5% (84.5% - 100%)
TRA (200)	Viewer C	+	10	20	8	0	0	95% (79.5% - 100%)
	Viewer B	-	0	0	1	19	21	100% (84.5% - 100%)
	Viewer A	+	0	0	2	19	21	100% (84.5% - 100%)
TRA (1000)	Viewer C	+	10	20	8	0	0	95% (79.5% - 100%)
	Viewer B	-	0	0	1	19	20	97.5% (84.5% - 100%)
	Viewer A	+	0	0	1	18	20	95% (84.5% - 100%)

Precision and Sensitivity

To investigate the precision and sensitivity each drug sample was analyzed at the following concentrations: cutoff +100%, cutoff -75%, cutoff -50%, cutoff -25%, cutoff -10%, cutoff -50%, cutoff +75% and the cutoff +100%. All concentrations were confirmed with GC-MS. The study was performed 2 runs /day and lasted 25 days using three different lots of the corresponding drug test. Totally 3 operators participated in the study of the corresponding drug test. Each of the 3 operators tests 2 aliquots at each concentration for each lot per day (2 runs/day), for a total of 50 determinations per concentration per lot of the corresponding drug test.

Drug Test	Approximate Concentration Sample (ng/mL)	Number of Determinations per Lot	Results Negative/Positive		
			Lot 1	Lot 2	Lot 3
			Lot 1	Lot 2	Lot 3
6-MAM	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	47/3	48/2	47/3
	12.5	50	5/45	6/44	6/44
	15	50	3/47	2/48	2/48
	17.5	50	0/50	0/50	0/50
	20	50	0/50	0/50	0/50
	22.5	50	0/50	0/50	0/50
	25	50	0/50	0/50	0/50
AMP (300)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
AMP (500)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
AMP (1000)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
EDDP (300)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
EDDP (500)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
AMP (1000)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
BAR	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
BUP (5)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
BUP (10)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
BZO (100)	2.5	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
	7.5	50	50/0	50/0	50/0
	10	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	17.5	50	50/0	50/0	50/0
	20	50	50/0	50/0	50/0
	22.5	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0

BZO (200)	175	50	0/50	0/50	0/50
	200	50	0/50	0/50	0/50
	225	50	0/50	0/50	0/50
	250	50	0/50	0/50	0/50
	275	50	0/50	0/50	0/50
	300	50	0/50	0/50	0/50
	325	50	0/50	0/50	0/50
	350	50	0/50	0/50	0/50
	375	50	0/50	0/50	0/50
	400	50	0/50	0/50	0/50
BZO (300)	0	50	50/0	50/0	50/0
	75	50	50/0	50/0	50/0
	150	50	50/0	50/0	50/0
	225	50	50/0	50/0	50/0
	300	50	6/44	5/45	6/44
	375	50	0/50	0/50	0/50
	450	50	0/50	0/50	0/50
	525	50	0/50	0/50	0/50
	600	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
COC (100)	25	50	50/0	50/0	50/0
	50	50	50/0	50/0	50/0
	75	50	50/0	50/0	50/0
	100	50	4/46	4/46	3/47
	125	50	0/50	0/50	0/50
	150	50	0/50	0/50	0/50
	175	50	0/50	0/50	0/50
	200	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	375	50	50/0	50/0	50/0
COC (150)	75	50	50/0	50/0	50/0
	112.5	50	50/0	50/0	50/0
	150	50	7/43	6/44	7/43
	187.5	50	0/50	0/50	0/50
	225	50	0/50	0/50	0/50
	262.5	50	0/50	0/50	0/50
	300	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	75	50	50/0	50/0	50/0
	150	50	50/0	50/0	50/0
COC (300)	225	50	50/0	50/0	50/0
	300	50	6/44	5/45	5/45
	375	50	0/50	0/50	0/50
	450	50	0/50	0/50	0/50
	525	50	0/50	0/50	0/50
	600	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	50	50	50/0	50/0	50/0
	100	50	50/0	50/0	50/0
	150	50	48/2	49/1	47/3
COT (200)	200	50	6/44	4/46	5/45
	250	50	4/46	3/47	2/48
	300	50	0/50	0/50	0/50
	350	50	0/50	0/50	0/50
	400	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	25	50	50/0	50/0	50/0
	50	50	50/0	50/0	50/0
	75	50	48/2	46/4	47/3
	100	50	6/44	5/45	5/45
EDDP (100)	125	50	2/48	3/47	5/45
	150	50	0/50	0/50	0/50
	175	50	0/50	0/50	0/50
	200	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	75	50	50/0	50/0	50/0
	150	50	50/0	50/0	50/0
	225	50	50/0	50/0	50/0
	300	50	6/44	5/45	6/44
	375	50	0/50	0/50	0/50
EDDP (300)	450	50	0/50	0/50	0/50
	525	50	0/50	0/50	0/50
	600	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	75	50	50/0	50/0	50/0
	150	50	50/0	50/0	50/0
	225	50	50/0	50/0	50/0
	300	50	50/0	50/0	50/0
	375	50	50/0	50/0	50/0
	450	50	50/0	50/0	50/0
ETG (300)	525	50	0/50	0/50	0/50
	600	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	125	50	50/0	50/0	50/0
	250	50	50/0	50/0	50/0
	375	50	50/0	50/0	50/0
	450	50	5/45	4/46	5/45
	525	50	0/50	0/50	0/50
	600	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
ETG (500)	125	50	50/0	50/0	50/0
	250	50	50/0	50/0	50/0
	375	50	50/0	50/0	50/0
	500	50	5/45	4/46	5/45
	625	50	0/50	0/50	0/50
	750	50	0/50	0/50	0/50
	875	50	0/50	0/50	0/50
	1000	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	5	50	50/0	50/0	50/0
FTY	10	50	50/0	50/0	50/0
	15	50	50/0	50/0	50/0
	20	50	4/46	5/45	5/45
	25	50	0/50	0/50	0/50
	30	50	0/50	0/50	0/50
	35	50	0/50	0/50	0/50
GAB	40	50	0/50	0/50	0/50
	0	50	50/0	50/0	50/0
	500	50	50/0	50/0	50/0
	1000	50	50/0	50/0	50/0
	1500	50	42/8	41/9	40/6
	2000	50	24/24	23/27	28/22
	2500	50	3/47	2/48	4/46
	3000	50	0/50	0/50	0/50
	3500	50	0/50	0/50	0/50