













PEN TEST-COCA-TEST **Please read instruction carefully before use.**

	 <p>Take all content out of the pouch.</p>		Hold the cap with your finger and shake well for 10 seconds.				
	Take the collection paper. Peel away the protective silicon cover. Hold it at the edge. Do not touch the exposed surface to avoid pre-contamination.		Remove cap, pointing the tube away from your body. Point the tube down and apply few drops on the paper to cover sampling area.				
	Wipe the suspected area, or touch the suspected substance with the test paper.		Observe for an immediate appearance of turquoise stains, indicating the presence of Cocaine, Crack and PCP. Check and compare results carefully with the color chart. Discard the test ingredients.				
	Holding the tube upright, squeeze the plastic tube to break the inside ampoule. Do not squeeze tube again over broken glass.		<table><tr><td>Cocaine</td><td rowspan="3">2-F </td></tr><tr><td>PCP</td></tr><tr><td>Crack</td></tr></table>	Cocaine	2-F 	PCP	Crack
Cocaine	2-F 						
PCP							
Crack							

Only an immediate color reaction should be considered as positive. Only colors similar to the color chart should be considered as positive results.
 For each test a new collection paper must be used. Use the collection paper enclosed in each pouch.
Caution: keep out of reach of children. In case of eye or skin contact, flush well with water.

1. This pouch contains one single test. Once ampoule is broken, the test should be used immediately.
2. Coca- Test detects the presence of: Cocaine, PCP and crack.
3. Coca-Test should only be used as a stand-alone test if one suspects one of these substances.
4. If one is not certain what the substance may be, or if there is a need to distinguish between cocaine/crack and PCP, one should use MD-1 for PCP and MD-2 to identify or detect cocaine or crack.
5. Use color chart to interpret the results.
6. Any positive color reactions may take 10 seconds or slightly longer to develop fully depending on the quantity and purity of the substance being tested.