











PEN TEST-HEROSOL

Please read instruction carefully before use.

 <p>①</p>	 <p>Take all content out of the pouch.</p>	 <p>⑤</p>	<p>Hold the cap with your finger and shake well for 10 seconds.</p>
 <p>②</p>	<p>Take the (#533) collection paper. Peel away the protective silicon cover. Hold it at the edge. Do not touch the exposed surface to avoid pre-contamination.</p>	 <p>⑥</p>	<p>Remove cap, pointing the tube away from your body. Point the tube down and apply few drops on the paper to cover sampling area.</p>
 <p>③</p>	<p>Wipe the suspected area, or touch the suspected substance with the test paper.</p>	 <p>⑦</p>	<p>Observe for an immediate appearance of purple stains, indicating the presence of Heroin, Methadone, Amphetamines, LSD and Methamphetamines. Check and compare results carefully with the color chart. Discard the test ingredients.</p>
 <p>④</p>	<p>Holding the tube upright, squeeze the plastic tube to break the inside ampoule. Do not squeeze tube again over broken glass.</p>		<p>Heroin LSD Methadone Methamphetamines Amphetamines</p> 

Only colors similar to the color chart should be considered as positive results.

For each test a new collection paper must be used. Use the collection paper enclosed in each pouch.

Caution: keep out of reach of children. In case of eye or skin contact, flush well with water.

1. This pouch contains two tests. Once ampoule is broken, the test should be used immediately.
2. Herosol detects the presence of: Heroin, Methadone, LSD and Methamphetamines.
3. If there is a need to distinguish between LSD, Heroin and Methamphetamines, one should use LSD test to identify or detect the presence of LSD, or use Meth-X to identify the presence of Methamphetamines.
4. Use color chart to interpret the results.
5. Any positive color reactions may take **10 seconds** or slightly longer to develop fully depending on the quantity and purity of the substance being tested.