














**PEN TEST-Opiatest**  
**Please read instruction carefully before use.**

①	 <p>Take all content out of the pouch.</p>	⑤	 <p>Hold the cap with your finger and shake well for 10 seconds.</p>
②	 <p>Take the collection paper. Peel away the protective silicon cover. Hold it at the edge. Do not touch the exposed surface to avoid pre-contamination.</p>	⑥	 <p>Remove cap, pointing the tube away from your body. Point the tube down and apply few drops on the paper to cover sampling area.</p>
③	 <p>Wipe the suspected area, or touch the suspected substance with the test paper.</p>	⑦	 <p>Check and compare results carefully with the color chart .</p> <p>Discard the test ingredients.</p>
④	 <p>Holding the tube upright, squeeze the plastic tube to break the inside ampoule. <b>Do not</b> squeeze tube again over broken glass.</p>		 <p>Opium      <b>2-B</b>      </p>
			 <p>Morphine      <b>2-A</b>      </p>
			 <p>Codeine      <b>1-E</b>      </p>

Only an immediate color reaction should be considered as positive. Unless the instructions for a particular drug type indicated differently. Only colors similar to the color chart should be considered as positive results. For each test a new collection paper must be used. Use the collection paper enclosed in each pouch.

**Caution:** keep out of reach of children. In case of exposure, flush with water.

1. This pouch contains one single test. Once ampoule is broken, the test should be used immediately.
2. This test identifies the presence of: Opium, Morphine and Codeine.
3. This test is used either with MD-1 as part of the narcotic screening test, or the test may be used as a stand-alone test if one suspects one of these substances.
4. Any positive color reactions may take 10 seconds or slightly longer to develop fully depending on the quantity and purity of the substance being tested.
5. A positive reaction may also look like a "circle" where the darkest part of the reaction may be on the perimeter of the color reaction.
6. Use color chart to interpret the results. Every reaction has a corresponding letter and number, for example: (2-A), Morphine.
7. Continue to the next step according to the color reaction received.
8. Use the next product as per the flow chart. The product name is indicated on the triangle sign.